

MTF1 antibody

Product Information

Catalog No.:	FNab05396
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Activates the metallothionein I promoter. Binds to the metal responsive element(MRE).

Immunogen information

Immunogen:	metal-regulatory transcription factor 1
Synonyms:	Metal regulatory transcription factor 1 MRE-binding transcription factor Transcription factor MTF-1 MTF1
Observed MW:	65-70 kDa
Uniprot ID :	Q14872

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab05396(MTF1 Antibody) at dilution of 1:1000