

## NUDT1 antibody

### Product Information

Catalog No.:	FNab05400
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	nudix(nucleoside diphosphate linked moiety X)-type motif 1
Synonyms:	Oxidized purine nucleoside triphosphate hydrolase 2-hydroxy-dATP diphosphatase 7,8-dihydro-8-oxoguanine triphosphatase 8-oxo-dGTPase Methylated purine nucleoside triphosphate hydrolase Nucleoside diphosphate-linked moiety X motif 1 (Nudix motif 1) NUDT1 MTH1
Observed MW:	18-24 kDa
Uniprot ID :	P36639

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab05400(NUDT1 antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab05400(NUDT1 antibody) at dilution of 1:300