

MTMR8 antibody

Product Information

Catalog No.:	FNab05415
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Phosphatase that acts on lipids with a phosphoinositol headgroup.

Immunogen information

Immunogen:	myotubularin related protein 8
Synonyms:	Myotubularin-related protein 8 Phosphatidylinositol-3,5-bisphosphate 3-phosphatase Phosphatidylinositol-3-phosphate phosphatase MTMR8
Observed MW:	55-58 kDa
Uniprot ID :	Q96EF0

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05415(MTMR8 antibody) at dilution of 1:300