

## MTX1 antibody

### Product Information

Catalog No.:	FNab05424
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


### Background

#### Immunogen information


Immunogen:	metaxin 1
Synonyms:	Metaxin-1 Mitochondrial outer membrane import complex protein 1 MTX1 MTX MTXN
Observed MW:	35 kDa
Uniprot ID :	Q13505

#### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon using FNab05424(MTX1 antibody) at dilution of 1:100



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab05424(MTX1 antibody) at dilution of 1:1000