

## PWWP3A antibody

### Product Information

Catalog No.:	FNab05437
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in the DNA damage response pathway by contributing to the maintenance of chromatin architecture. Recruited to the vicinity of DNA breaks by TP53BP1 and plays an accessory role to facilitate damage-induced chromatin changes and promoting chromatin relaxation. Required for efficient DNA repair and cell survival following DNA damage.

### Immunogen information


Immunogen:	melanoma associated antigen(mutated) 1
Synonyms:	PWWP domain-containing DNA repair factor 3A (PWWP3A) Mutated melanoma-associated antigen 1 (MUM-1) PWWP domain-containing protein MUM1 Protein expandere PWWP3A EXPAND1 MUM1
Observed MW:	70-80 kDa
Uniprot ID :	Q2TAK8

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human endometrial cancer using FNab05437(MUM1 antibody) at dilution of 1:50



HT-1080 cells were subjected to SDS PAGE followed by western blot with FNab05437(MUM1 antibody) at dilution of 1:300