

CAVIN4 antibody

Product Information

Catalog No.:	FNab05440
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Induces RHOA activation and activates NPPA transcription and myofibrillar organization through the Rho/ROCK signaling pathway.

Immunogen information

Immunogen:	muscle-related coiled-coil protein
Synonyms:	Caveolae-associated protein 4 Muscle-related coiled-coil protein Muscle-restricted coiled-coil protein CAVIN4 MURC
Observed MW:	48-53 kDa
Uniprot ID :	Q5BKC8

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab05440(MURC Antibody) at dilution of 1:600