

MKLN1 antibody

Product Information

| | |
|---------------|----------------------------------------------------------------------------------------------------------------|
| Catalog No.: | FNab05443 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


None

Immunogen information


| | |
|--------------|------------------------------------------------------------|
| Immunogen: | muskelin 1, intracellular mediator containing kelch motifs |
| Synonyms: | Muskelin MKLN1 |
| Observed MW: | 72 kDa |
| Uniprot ID : | Q9UL63 |

Application

| | |
|-----------------------|-----------------------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB, IHC |
| Recommended dilution: | WB: 1:200-1:1000; IHC: 1:20-1:200 |
| Image: | |



Immunohistochemistry of paraffin-embedded human gliomas using FNab05443(MKLN1 antibody) at dilution of 1:50



Y79 cells were subjected to SDS PAGE followed by western blot with FNab05443(MKLN1 antibody) at dilution of 1:1000