

MYBBP1A antibody

Product Information

Catalog No.:	FNab05453
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


This gene encodes a nucleolar transcriptional regulator that was first identified by its ability to bind specifically to the Myb proto-oncogene protein. The encoded protein is thought to play a role in many cellular processes including response to nucleolar stress, tumor suppression and synthesis of ribosomal DNA. Alternate splicing results in multiple transcript variants.

Immunogen information


Immunogen:	MYB binding protein (P160) 1a
Synonyms:	Myb-binding protein 1A MYBBP1A P160
Observed MW:	150 kDa
Uniprot ID :	Q9BQG0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human stomach using FNab05453(MYBBP1A antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05453(MYBBP1A antibody) at dilution of 1:1000