

MYCBPAP antibody

Product Information

Catalog No.:	FNab05458
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in spermatogenesis. May be involved in synaptic processes(By similarity).

Immunogen information

Immunogen:	MYCBP associated protein
Synonyms:	MYCBP-associated protein AMAM-1 AMY-1-binding protein 1 (AMAP-1) MYCBPAP AMAP1
Observed MW:	108 kDa
Uniprot ID :	Q8TBZ2

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab05458(MYCBPAP antibody) at dilution of 1:300