

MYCL antibody

Product Information

Catalog No.:	FNab05459
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	v-myc myelocytomatosis viral oncogene homolog 1, lung carcinoma derived(avian)
Synonyms:	Protein L-Myc Class E basic helix-loop-helix protein 38 (bHLHe38) Protein L-Myc-1 V-myc myelocytomatosis viral oncogene homolog MYCL BHLHE38 LMYC MYCL1
Observed MW:	35 kDa
Uniprot ID :	P12524

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:1000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05459(MYCL1 antibody) at dilution of 1:500