

MYEOV antibody

Product Information

Catalog No.:	FNab05466
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	myeloma overexpressed(in a subset of t(11;14) positive multiple myelomas)
Synonyms:	Myeloma-overexpressed gene protein Oncogene in multiple myeloma MYEOV OCIM
Observed MW:	55-60 kDa
Uniprot ID :	Q96EZ4

Application

Reactivity:	Human, rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab05466(MYEOV antibody) at dilution of 1:50



SGC-7901 cells were subjected to SDS PAGE followed by western blot with FNab05466(MYEOV antibody) at dilution of 1:400