

MYH3 antibody

Product Information

Catalog No.:	FNab05475
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Muscle contraction.

Immunogen information

Immunogen:	myosin, heavy chain 3, skeletal muscle, embryonic
Synonyms:	Myosin-3 Muscle embryonic myosin heavy chain Myosin heavy chain 3 Myosin heavy chain, fast skeletal muscle, embryonic SMHCE MYH3
Observed MW:	220-230 kDa
Uniprot ID :	P11055

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse embryo tissue were subjected to SDS PAGE followed by western blot with FNab05475(MYH3 Antibody) at dilution of 1:1000