

MYH4-Specific antibody

Product Information

Catalog No.:	FNab05476
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	myosin, heavy chain 4, skeletal muscle
Synonyms:	Myosin-4 Myosin heavy chain 2b (MyHC-2b) Myosin heavy chain 4 Myosin heavy chain IIb (MyHC-IIb) Myosin heavy chain, skeletal muscle, fetal MYH4
Observed MW:	230 kDa
Uniprot ID :	Q9Y623

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using FNab05476(MYH4-Specific Antibody) at dilution of 1:50



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab05476(MYH4-Specific Antibody) at dilution of 1:500