

MYL4 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab05488 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background

Regulatory light chain of myosin. Does not bind calcium.

Immunogen information

| | |
|--------------|---|
| Immunogen: | myosin, light chain 4, alkali; atrial, embryonic |
| Synonyms: | Myosin light chain 4 Myosin light chain 1, embryonic muscle/atrial isoform Myosin light chain alkali GT-1 isoform MYL4 MLC1 |
| Observed MW: | 22 kDa |
| Uniprot ID : | P12829 |

Application

| | |
|-----------------------|-----------------------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB, IHC |
| Recommended dilution: | WB: 1:500-1:2000; IHC: 1:20-1:200 |
| Image: | |



human liver tissue were subjected to SDS PAGE followed by western blot with FNab05488(MYL4 antibody) at dilution of 1:1000