

## MYLK2 antibody

### Product Information

Catalog No.:	FNab05494
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Implicated in the level of global muscle contraction and cardiac function. Phosphorylates a specific serine in the N-terminus of a myosin light chain.

### Immunogen information

Immunogen:	myosin light chain kinase 2
Synonyms:	Myosin light chain kinase 2, skeletal/cardiac muscle (MLCK2) MYLK2
Observed MW:	65 kDa
Uniprot ID :	Q9H1R3

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab05494(MYLK2 antibody) at dilution of 1:800