

## MYOM1-Specific antibody

### Product Information

Catalog No.:	FNab05516
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	myomesin 1, 185 kDa
Synonyms:	Myomesin-1 190 kDa connectin-associated protein 190 kDa titin-associated protein Myomesin family member 1 MYOM1
Observed MW:	190 kDa
Uniprot ID :	P52179

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using FNab05516(MYOM1-Specific Antibody) at dilution of 1:50



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab05516(MYOM1-Specific antibody) at dilution of 1:600