

SSNA1 antibody

Product Information

Catalog No.:	FNab05534
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	Sjogren syndrome nuclear autoantigen 1
Synonyms:	Microtubule nucleation factor SSNA1 Nuclear autoantigen of 14 kDa Sjogren syndrome nuclear autoantigen 1 SSNA1 NA14
Observed MW:	14 kDa
Uniprot ID :	O43805

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab05534(NA14 antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab05534(NA14 antibody) at dilution of 1:400