

NAP1L4 antibody

Product Information

Catalog No.:	FNab05545
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	nucleosome assembly protein 1-like 4
Synonyms:	Nucleosome assembly protein 1-like 4 Nucleosome assembly protein 2 (NAP-2) NAP1L4 NAP2
Observed MW:	52 kDa
Uniprot ID :	Q99733

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab05545(NAP1L4 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05545(NAP1L4 antibody) at dilution of 1:1000