

## NAPRT antibody

### Product Information

Catalog No.:	FNab05549
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	0C6
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	nicotinate phosphoribosyltransferase domain containing 1
Synonyms:	Nicotinate phosphoribosyltransferase (NAPRTase) FHA-HIT-interacting protein Nicotinate phosphoribosyltransferase domain-containing protein 1 NAPRT FHIP NAPRT1
Observed MW:	51 kDa
Uniprot ID :	Q6XQN6

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using FNab05549(NAPRT1 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05549(NAPRT1 Antibody) at dilution of 1:1500