

NAA50 antibody

Product Information

Catalog No.:	FNab05563
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	N-acetyltransferase 13(GCN5-related)
Synonyms:	N-alpha-acetyltransferase 50 (hNaa50p) N-acetyltransferase 13 N-acetyltransferase 5 (hNAT5) N-acetyltransferase san homolog (hSAN) N-epsilon-acetyltransferase 50 NatE catalytic subunit NAA50 MAK3 NAT13 NAT5
Observed MW:	19 kDa
Uniprot ID :	Q9GZZ1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab05563(NAT13 antibody) at dilution of 1:300