

NBAS antibody

Product Information

Catalog No.:	FNab05565
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in Golgi-to-endoplasmic reticulum(ER) retrograde transport; the function is proposed to depend on its association in the NRZ complex which is believed to play a role in SNARE assembly at the ER(PubMed:19369418).

Immunogen information

Immunogen:	neuroblastoma amplified sequence
Synonyms:	NBAS subunit of NRZ tethering complex Neuroblastoma-amplified gene protein Neuroblastoma-amplified sequence NBAS NAG
Observed MW:	250 kDa
Uniprot ID :	A2RRP1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab05565(NABS antibody) at dilution of 1:500