

NCALD antibody

Product Information

Catalog No.:	FNab05572
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in the calcium-dependent regulation of rhodopsin phosphorylation. Binds three calcium ions.

Immunogen information

Immunogen:	neurocalcin delta
Synonyms:	Neurocalcin-delta NCALD
Observed MW:	18-22 kDa
Uniprot ID :	P61601

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab05572(NCALD antibody) at dilution of 1:2000