

NDUFAF2 antibody

Product Information

Catalog No.:	FNab05615
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Acts as a molecular chaperone for mitochondrial complex I assembly.

Immunogen information


Immunogen:	NADH dehydrogenase(ubiquinone) 1 alpha subcomplex, assembly factor 2
Synonyms:	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 2 B17.2-like (B17.2L) Mimitin Myc-induced mitochondrial protein (MMTN) NDUFA12-like protein NDUFAF2 NDUFA12L
Observed MW:	25 kDa
Uniprot ID :	Q8N183

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab05615(NDUF2 antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab05615(NDUF2 antibody) at dilution of 1:800