

NDUFAF4 antibody

Product Information

Catalog No.:	FNab05616
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the assembly of mitochondrial NADH:ubiquinone oxidoreductase complex(complex I). May be involved in cell proliferation and survival of hormone-dependent tumor cells. May be a regulator of breast tumor cell invasion.

Immunogen information

Immunogen:	NADH dehydrogenase(ubiquinone) 1 alpha subcomplex, assembly factor 4
Synonyms:	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 4 Hormone-regulated proliferation-associated protein of 20 kDa NDUFAF4 C6orf66 HRPAP20
Observed MW:	21 kDa
Uniprot ID :	Q9P032

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab05616(NDUFAF4 Antibody) at dilution of 1:1200