

## NDUFB10 antibody

### Product Information

Catalog No.:	FNab05618
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


### Background

#### Immunogen information


Immunogen:	NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 10, 22kDa
Synonyms:	NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 10 Complex I-PDSW (CI-PDSW) NADH-ubiquinone oxidoreductase PDSW subunit NDUFB10
Observed MW:	21 kDa
Uniprot ID :	O96000

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded mouse stomach using FNab05618(NDUF B10 antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab05618(NDUF B10 antibody) at dilution of 1:1000