

NSMF antibody

Product Information

Catalog No.:	FNab05655
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	nasal embryonic LHRH factor
Synonyms:	NMDA receptor synaptonuclear signaling and neuronal migration factor Nasal embryonic luteinizing hormone-releasing hormone factor (Nasal embryonic LHRH factor) NSMF NELF
Observed MW:	54-60 kDa
Uniprot ID :	Q6X4W1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab05655(NELF antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05655(NELF Antibody) at dilution of 1:1000