

NELFA antibody

Product Information


Catalog No.:	FNab05656
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	Wolf-Hirschhorn syndrome candidate 2
Synonyms:	Negative elongation factor A (NELF-A) Wolf-Hirschhorn syndrome candidate 2 protein NELFA WHSC2
Observed MW:	57 kDa
Uniprot ID :	Q9H3P2

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IF, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:1000; IF: 1:10-1:100
Image:	



IP Result of anti-NELF-A (IP:FNab05656, 4ug;
Detection:FNab05656 1:500) with mouse brain
tissue lysate 3440ug.



HeLa cells were subjected to SDS PAGE followed
by western blot with FNab05656(NELF-A
Antibody) at dilution of 1:1000