

NELL1 antibody

Product Information

Catalog No.:	FNab05659
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a role in the control of cell growth and differentiation. Promotes osteoblast cell differentiation and terminal mineralization.

Immunogen information

Immunogen:	NEL-like 1
Synonyms:	Protein kinase C-binding protein NELL1 NEL-like protein 1 Nel-related protein 1 NELL1 NRP1
Observed MW:	84 kDa
Uniprot ID :	Q92832

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse pancreas tissue were subjected to SDS
PAGE followed by western blot with
FNab05659(NELL1 antibody) at dilution of 1:500