

## NEO1-Specific antibody

### Product Information

Catalog No.:	FNab05662
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	neogenin
Synonyms:	Neogenin Immunoglobulin superfamily DCC subclass member 2 NEO1 IGDCC2 NGN
Observed MW:	160 kDa
Uniprot ID :	Q92859

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab05662(Neogenin-Specific antibody) at dilution of 1:500