

NEUROD4 antibody

Product Information

Catalog No.:	FNab05675
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Probably acts as a transcriptional activator. Mediates neuronal differentiation. Required for the regulation of amacrine cell fate specification in the retina(By similarity).

Immunogen information

Immunogen:	neurogenic differentiation 4
Synonyms:	Neurogenic differentiation factor 4 (NeuroD4) Class A basic helix-loop-helix protein 4 (bHLHa4) Protein atonal homolog 3 (ATH-3, Atoh3) NEUROD4 ATH3 ATOH3 BHLHA4
Observed MW:	50-55 kDa
Uniprot ID :	Q9HD90

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Y79 cells were subjected to SDS PAGE followed by western blot with FNab05675(NEUROD4 Antibody) at dilution of 1:1000