

## ILF2 antibody

### Product Information

Catalog No.:	FNab05687
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	interleukin enhancer binding factor 2, 45kDa
Synonyms:	Interleukin enhancer-binding factor 2 Nuclear factor of activated T-cells 45 kDa ILF2 NF45
Observed MW:	43 kDa
Uniprot ID :	Q12905

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer slide using FNab05687(NF45 antibody) at dilution of 1:100



K-562 cells were subjected to SDS PAGE followed by western blot with FNab05687(NF45 antibody) at dilution of 1:500