

NFKB1 p105,p50-Specific antibody

Product Information

Catalog No.:	FNab05704
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	nuclear factor of kappa light polypeptide gene enhancer in B-cells 1
Synonyms:	Nuclear factor NF-kappa-B p105 subunit DNA-binding factor KBF1 EBP-1 Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1 Nuclear factor NF-kappa-B p50 subunit NFKB1
Observed MW:	50 kDa, 105 kDa
Uniprot ID :	P19838

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



A431 cells were subjected to SDS PAGE followed by western blot with FNab05704(NFKB1,p105,p50-Specific antibody) at dilution of 1:200