

## NFKBIE antibody

### Product Information

Catalog No.:	FNab05708
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Inhibits NF-kappa-B by complexing with and trapping it in the cytoplasm. Inhibits DNA-binding of NF-kappa-B p50-p65 and p50-c-Rel complexes.

### Immunogen information

Immunogen:	nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor, epsilon
Synonyms:	NF-kappa-B inhibitor epsilon (NF-kappa-BIE) I-kappa-B-epsilon (IkB-E, IkB-epsilon, IkappaBepsilon) NFKBIE IKBE
Observed MW:	51 kDa
Uniprot ID :	O00221

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



human spleen tissue were subjected to SDS PAGE followed by western blot with FNab05708(NFKBIE antibody) at dilution of 1:500