

NEFM-Specific antibody

Product Information

Catalog No.:	FNab05712
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	neurofilament, medium polypeptide
Synonyms:	Neurofilament medium polypeptide (NF-M) 160 kDa neurofilament protein Neurofilament 3 Neurofilament triplet M protein NEFM NEF3 NFM
Observed MW:	150 kDa
Uniprot ID :	P07197

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, WB, IHC
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:50-1:500; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab05712 (NEFM-Specific Antibody) at dilution of 1:200



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab05712 (NF-M-Specific antibody) at dilution of 1:1000