

NFYC antibody

Product Information

Catalog No.:	FNab05717
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Component of the sequence-specific heterotrimeric transcription factor(NF-Y) which specifically recognizes a 5'-CCAAT-3' box motif found in the promoters of its target genes. NF-Y can function as both an activator and a repressor, depending on its interacting cofactors.

Immunogen information

Immunogen:	nuclear transcription factor Y, gamma
Synonyms:	Nuclear transcription factor Y subunit gamma CAAT box DNA-binding protein subunit C Nuclear transcription factor Y subunit C (NF-YC) Transactivator HSM-1/2 NFYC
Observed MW:	48 kDa
Uniprot ID :	Q13952

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab05717(NFYC antibody) at dilution of 1:1000