

NIF3L1 antibody

Product Information

Catalog No.:	FNab05732
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May function as a transcriptional corepressor through its interaction with COPS2, negatively regulating the expression of genes involved in neuronal differentiation.

Immunogen information

Immunogen:	NIF3 NGG1 interacting factor 3-like 1(S. pombe)
Synonyms:	NIF3-like protein 1 Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 1 protein NIF3L1 ALS2CR1
Observed MW:	45 kDa
Uniprot ID :	Q9GZT8

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab05732(NIF3L1 Antibody) at dilution of 1:600