

NINJ2 antibody

Product Information

Catalog No.:	FNab05734
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	7A2
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Homophilic cell adhesion molecule that promotes axonal growth. May play a role in nerve regeneration and in the formation and function of other tissues.

Immunogen information

Immunogen:	ninjurin 2
Synonyms:	Ninjurin-2 Nerve injury-induced protein 2 NINJ2
Observed MW:	None
Uniprot ID :	Q9NZG7

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab05734(NINJ2 Antibody) at dilution of 1:50