

## PSME3IP1 antibody

### Product Information

Catalog No.:	FNab05735
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	NEFA-interacting nuclear protein NIP30
Synonyms:	PSME3-interacting protein NEFA-interacting nuclear protein NIP30 PA28G-interacting protein PSME3IP1 C16orf94 FAM192A NIP30 PIP30
Observed MW:	29 kDa
Uniprot ID :	Q9GZU8

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05735(NIP30 antibody) at dilution of 1:500