

KLRC1 antibody

Product Information

Catalog No.:	FNab05744
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	killer cell lectin-like receptor subfamily C, member 1
Synonyms:	NKG2-A/NKG2-B type II integral membrane protein CD159 antigen-like family member A NK cell receptor A NKG2-A/B-activating NK receptor KLRC1 NKG2A
Observed MW:	33 kDa
Uniprot ID :	P26715

Application

Reactivity:	Human
Tested Application:	ELISA,WB
Recommended dilution:WB:	1:500-1:1000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab05744(NKG2A antibody) at dilution of 1:800