

NKX1-2 antibody

Product Information

Catalog No.:	FNab05749
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May function in cell specification, particularly in the CNS.

Immunogen information


Immunogen:	NK1 homeobox 2
Synonyms:	NK1 transcription factor-related protein 2 Homeobox protein SAX-1 NKX-1.1 NKX1-2 C10orf121 NKX1.1
Observed MW:	32 kDa
Uniprot ID :	Q9UD57

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab05749(NKX1-2 Antibody) at dilution of 1:50



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab05749(NKX1-2 Antibody) at dilution of 1:300