

## NMD3 antibody

### Product Information

Catalog No.:	FNab05760
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Acts as an adapter for the XPO1/CRM1-mediated export of the 60S ribosomal subunit.

### Immunogen information


Immunogen:	NMD3 homolog
Synonyms:	60S ribosomal export protein NMD3 (hNMD3) NMD3
Observed MW:	58 kDa
Uniprot ID :	Q96D46

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-NMD3 (IP:FNab05760, 3ug;  
Detection:FNab05760 1:800) with HepG2 cells  
lysate 1720ug.



HepG2 cells were subjected to SDS PAGE  
followed by western blot with FNab05760(NMD3  
antibody) at dilution of 1:1000