

RIOX1 antibody

Product Information

Catalog No.:	FNab05776
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	chromosome 14 open reading frame 169
Synonyms:	Ribosomal oxygenase 1 60S ribosomal protein L8 histidine hydroxylase Bifunctional lysine-specific demethylase and histidyl-hydroxylase NO66 Myc-associated protein with JmjC domain Nucleolar protein 66 (hsNO66) Ribosomal oxygenase NO66 (ROX) RIOX1 C14orf169 MAPJD NO66
Observed MW:	71 kDa
Uniprot ID :	Q9H6W3

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



A431 cells were subjected to SDS PAGE followed by western blot with FNab05776(NO66/C14orf169 antibody) at dilution of 1:600