

NOB1 antibody

Product Information

Catalog No.:	FNab05778
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in mRNA degradation.

Immunogen information

Immunogen:	NIN1/RPN12 binding protein 1 homolog
Synonyms:	RNA-binding protein NOB1 Phosphorylation regulatory protein HP-10 Protein ART-4 NOB1 ART4 NOB1P PSMD8BP1
Observed MW:	47 kDa
Uniprot ID :	Q9ULX3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab05778(NOB1 antibody) at dilution of 1:500