

NOB1 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab05779 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Protein A+G purification |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Mouse |
| Clonality: | monoclonal |
| Clone ID: | 8G0 |
| IsoType: | IgG2a |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


May play a role in mRNA degradation.

Immunogen information


| | |
|--------------|--|
| Immunogen: | NIN1/RPN12 binding protein 1 homolog |
| Synonyms: | RNA-binding protein NOB1 Phosphorylation regulatory protein HP-10 Protein ART-4 NOB1 ART4 NOB1P PSMD8BP1 |
| Observed MW: | 47 kDa |
| Uniprot ID : | Q9ULX3 |

Application

| | |
|-----------------------|------------------------------------|
| Reactivity: | Human |
| Tested Application: | ELISA, WB, IHC |
| Recommended dilution: | WB: 1:1000-1:6000; IHC: 1:20-1:200 |
| Image: | |



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab05779(NOB1 Antibody) at dilution of 1:200



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab05779(NOB1 Antibody) at dilution of 1:3000