

NOC4L antibody

Product Information

Catalog No.:	FNab05780
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	nucleolar complex associated 4 homolog
Synonyms:	Nucleolar complex protein 4 homolog (NOC4 protein homolog) NOC4-like protein Nucleolar complex-associated protein 4-like protein NOC4L
Observed MW:	59 kDa
Uniprot ID :	Q9BVI4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05780(NOC4L antibody) at dilution of 1:300