

NOL8 antibody

Product Information

Catalog No.:	FNab05788
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays an essential role in the survival of diffuse-type gastric cancer cells. Acts as a nucleolar anchoring protein for DDX47. May be involved in regulation of gene expression at the post-transcriptional level or in ribosome biogenesis in cancer cells.

Immunogen information

Immunogen:	nucleolar protein 8
Synonyms:	Nucleolar protein 8 Nucleolar protein Nop132 NOL8 C9orf34 NOP132
Observed MW:	150 kDa
Uniprot ID :	Q76FK4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



SGC-7901 cells were subjected to SDS PAGE followed by western blot with FNab05788(NOL8 antibody) at dilution of 1:500