

NOP56 antibody

Product Information

Catalog No.:	FNab05795
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the early to middle stages of 60S ribosomal subunit biogenesis. Core component of box C/D small nucleolar ribonucleoprotein(snoRNP) particles. Required for the biogenesis of box C/D snoRNAs such U3, U8 and U14 snoRNAs.

Immunogen information

Immunogen:	NOP56 ribonucleoprotein homolog(yeast)
Synonyms:	Nucleolar protein 56 Nucleolar protein 5A NOP56 NOL5A
Observed MW:	58-66 kDa
Uniprot ID :	O00567

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab05795(NOP56 antibody) at dilution of 1:1000