

CCN3 antibody

Product Information

Catalog No.:	FNab05801
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Immediate-early protein playing a role in various cellular processes including proliferation, adhesion, migration, differentiation and survival(PubMed:15181016, PubMed:15611078, PubMed:12695522, PubMed:21344378, PubMed:12050162). Acts by binding to integrins or membrane receptors such as NOTCH1(PubMed:12695522, PubMed:21344378, PubMed:15611078). Essential regulator of hematopoietic stem and progenitor cell function(PubMed:17463287). Inhibits myogenic differentiation through the activation of Notch-signaling pathway(PubMed:12050162). Inhibits vascular smooth muscle cells proliferation by increasing expression of cell-cycle regulators such as CDKN2B or CDKN1A independently of TGFB1 signaling(PubMed:20139355). Ligand of integrins ITGAV:ITGB3 and ITGA5:ITGB1, acts directly upon endothelial cells to stimulate pro-angiogenic activities and induces angiogenesis. In endothelial cells, supports cell adhesion, induces directed cell migration(chemotaxis) and promotes cell survival(PubMed:12695522). Plays also a role in cutaneous wound healing acting as integrin receptor ligand. Supports skin fibroblast adhesion through ITGA5:ITGB1 and ITGA6:ITGB1 and induces fibroblast chemotaxis through ITGAV:ITGB5. Seems to enhance bFGF-induced DNA synthesis in fibroblasts(PubMed:15611078). Involved in bone regeneration as a negative regulator(By similarity). Enhances the articular chondrocytic phenotype, whereas it repressed the one representing endochondral ossification(PubMed:21871891). Impairs pancreatic beta-cell function, inhibits beta-cell proliferation and insulin secretion(By similarity). Plays a role as negative regulator of endothelial pro-inflammatory activation reducing monocyte adhesion, its anti-inflammatory effects occur secondary to the inhibition of NF-kappaB signaling pathway(PubMed:21063504). Contributes to the control and coordination of inflammatory processes in atherosclerosis(By similarity). Attenuates inflammatory pain through regulation of IL1B-and TNF-induced MMP9, MMP2 and CCL2 expression. Inhibits MMP9 expression through ITGB1 engagement(PubMed:21871891).

Immunogen information


Immunogen: nephroblastoma overexpressed gene
 Synonyms: CCN family member 3|Cellular communication network factor 3|Insulin-like growth factor-binding protein 9 (IBP-9, IGF-binding protein 9, IGFBP-9)|Nephro blastoma-overexpressed gene protein homolog|Protein NOV homolog (NovH)|CCN3|IGFBP9|NOV|NOVH
 Observed MW: 45 kDa
 Uniprot ID : P48745

Application

Reactivity: Human, Mouse, Rat
 Tested Application: ELISA, WB, IHC
 Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200
 Image:



Immunohistochemistry of paraffin-embedded human nephroblastoma using FNab05801(NOV antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab05801(NOV antibody) at dilution of 1:500