

NOVA2 antibody

Product Information

Catalog No.:	FNab05802
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May regulate RNA splicing or metabolism in a specific subset of developing neurons(By similarity). Binds single strand RNA.

Immunogen information

Immunogen:	neuro-oncological ventral antigen 2
Synonyms:	RNA-binding protein Nova-2 Astrocytic NOVA1-like RNA-binding protein Neuro-oncological ventral antigen 2 NOVA2 ANOVA NOVA3
Observed MW:	60-65 kDa
Uniprot ID :	Q9UNW9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab05802(NOVA2 antibody) at dilution of 1:500