

## NOX3-Specific antibody

### Product Information

Catalog No.:	FNab05805
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	NADPH oxidase 3
Synonyms:	NADPH oxidase 3 Mitogenic oxidase 2 (MOX-2) gp91phox homolog 3 (GP91-3) NOX3 MOX2
Observed MW:	52 kDa
Uniprot ID :	Q9HBY0

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab05805(NOX3-Specific antibody) at dilution of 1:500